

In the Specification:

Please amend the specification as follows:

Paragraph bridging pages 6 and 7:

The pre-mounted implant device is delivered ~~steril~~ sterile in a plastic package 11 comprising a titanium packaging sleeve 12 in order to retain the implant device in a predetermined position in the plastic package, see figure 3. At the surgery the plastic package is broken by removing the plastic lid 17 and the pre-mounted implant device is then separated from the titanium packaging sleeve 12 by means of the tool 4 and ~~said~~ the lifting function.

Page 7, first paragraph:

By placing the ~~pre-mounted~~ pre-mounted implant device in a titanium packaging sleeve 12 it is protected so that the tool will not come into contact with the fixture part when the implant device is removed from the packaging. A sealing ring 18 is arranged on the cylindrical outer surface of the plastic package to provide a tightening between the plastic package and the lid 17. The sealing ring 18 can be adjusted in the longitudinal direction to provide a tightening even for different positions of the plastic lid 17 on the package.